

PJA2 Polyclonal Antibody

Catalog No.	A6444	Category	Polyclonal Antibodies
Applications	WB, IF	Observed MW	Refer to Figures
Cross-reactivity	Human, Mouse, Rat	Calculated MW	76kDa/78kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 439-708 of human PJA2 (NP_055634.3).
Gene ID	9867
Swiss prot	O43164
Synonyms	PJA2; Neurodap1; RNF131; E3 ubiquitin-protein ligase Praja-2

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000

IF 1:50 -
1:100



Immunofluorescence - PJA2 Polyclonal Antibody (A6444)