

SLC19A1 Polyclonal Antibody

Catalog No.	A12819	Category	Polyclonal Antibodies
Applications	WB	Observed MW	65kDa
Cross-reactivity	Human	Calculated MW	54kDa/60kDa/64kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 452-591 of human SLC19A1 (NP_919231.1).
Gene ID	6573
Swiss prot	P41440
Synonyms	SLC19A1; CHMD; FOLT; IFC1; REFC; RFC1; folate transporter 1

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

The membrane protein encoded by this gene is a transporter of folate and is involved in the regulation of intracellular concentrations of folate. Three transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 -
1:1000



Western blot - SLC19A1 Polyclonal Antibody (A12819)