

CD300C Polyclonal Antibody

Catalog No.	A10420	Category	Polyclonal Antibodies
Applications	WB	Observed MW	52kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	24kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-183 of human CD300C (NP_006669.1).
Gene ID	10871
Swiss prot	Q08708
Synonyms	CD300C; CLM-6; CMRF-35; CMRF-35A; CMRF35; CMRF35-A1; CMRF35A; CMRF35A1; IGSF16; LIR; CD300c molecule


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - CD300C Polyclonal Antibody (A10420)