

REEP5 Polyclonal Antibody

Catalog No.	A10392	Category	Polyclonal Antibodies
Applications	WB	Observed MW	21kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	14kDa/21kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-189 of human REEP5 (NP_005660.4).
Gene ID	7905
Swiss prot	Q00765
Synonyms	REEP5; C5orf18; D5S346; DP1; TB2; YOP1; Yip2e; receptor accessory protein 5


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - REEP5 Polyclonal Antibody (A10392)