

NADK Polyclonal Antibody

Catalog No.	A9111	Category	Polyclonal Antibodies
Applications	WB	Observed MW	49kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	46kDa/49kDa/63kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human NADK (NP_075394.3).
Gene ID	65220
Swiss prot	O95544
Synonyms	NADK; dj283E3.1; NAD kinase


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:1000 -
1:2000



Western blot - NADK Polyclonal Antibody (A9111)