

SLC8A1 Polyclonal Antibody

Catalog No.	A7419	Category	Polyclonal Antibodies
Applications	WB	Observed MW	109kDa
Cross-reactivity	Human	Calculated MW	104kDa/107kDa/108kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-520 of human SLC8A1 (NP_066920.1).
Gene ID	6546
Swiss prot	P32418
Synonyms	SLC8A1; NCX1; sodium/calcium exchanger 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - SLC8A1 Polyclonal Antibody (A7419)