

## Anti-CD38 Antibody [ABT-CD38] (A100030)

### Specifications:

Name:	Anti-CD38 Antibody [ABT-CD38]
Description:	Mouse monoclonal (ABT-CD38) antibody to CD38.
Applications:	IHC
Reactivity:	Human
Immunogen:	Synthetic peptide derived from human CD38.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	ABT-CD38
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Immunogen affinity purified.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline with 0.5% Bovine Serum Albumin and 0.1% Kathon CG/ICP.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.