

ZNF22 Polyclonal Antibody

Catalog No.	A14798	Category	Polyclonal Antibodies
Applications	WB	Observed MW	26kDa
Cross-reactivity	Mouse	Calculated MW	25kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-224 of human ZNF22 (NP_008894.2).
Gene ID	7570
Swiss prot	P17026
Synonyms	ZNF22; HKR-T1; KOX15; ZNF422; Zfp422; zinc finger protein 22


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - ZNF22 Polyclonal Antibody (A14798)