


## Anti-Apolipoprotein E Antibody (A84074)

### Specifications:

Name:	Anti-Apolipoprotein E Antibody
Description:	Goat polyclonal antibody to Apolipoprotein E.
Applications:	ELISA, WB
Reactivity:	Human
Immunogen:	Synthetic peptide corresponding to Human APOE (C terminal).
Sequence:	C-VG TSAAPVPSDNH
Host:	Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Concentration:	100 µg at 0.5 mg/ml.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 7.30, with 0.02% Sodium Azide and 0.5% BSA.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

### Images:



Anti-APOE Antibody (A84074) (2.5µg/ml) staining of paraffin embedded Human Liver. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.

## Anti-Apolipoprotein E Antibody (A84074)

Images continued:

250kDa  
150kDa  
100kDa  
75kDa  
50kDa  
37kDa  
25kDa  
20kDa  
15kDa

Anti-APOE Antibody (A84074) (2 $\mu$ g/ml) staining of Human Brain (Hippocampus) lysate (35 $\mu$ g protein in RIPA buffer). Detected by chemiluminescence.