

LRRTM4 Polyclonal Antibody

Catalog No.	A15898	Category	Polyclonal Antibodies
Applications	WB	Observed MW	80kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	59kDa/67kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-420 of human LRRTM4 (NP_001269857.1).
Gene ID	80059
Swiss prot	Q86VH4
Synonyms	LRRTM4; leucine rich repeat transmembrane neuronal 4


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - LRRTM4 Polyclonal Antibody (A15898)