

NAIF1 Polyclonal Antibody

Catalog No.	A15217	Category	Polyclonal Antibodies
Applications	WB	Observed MW	35kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	17kDa/35kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human NAIF1 (NP_931045.1).
Gene ID	203245
Swiss prot	Q69YI7
Synonyms	NAIF1; C9orf90; bA379C10.2; nuclear apoptosis-inducing factor 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - NAIF1 Polyclonal Antibody (A15217)