

## Declaration of Conformity

For the: Variable Heat Sealer



The Variable Heat Sealer satisfies the following standards and directives:

- IEC 61010:2010 satisfying the Low Voltage Directive 2014/35/EU
- EN 61326:2013 satisfying the EMC directive 2004/108/EC
- EN 61010-2-010:2014 satisfying the Low Voltage Directive 2014/35/EU
- BSEN 50581:2012 satisfying the RoHS directive 2011/65/EU
- UL 61010-1:2001 2nd Edition/CAN C22.2 CSA 61010-1
- FCC CFR 47 Parts 15.107 and 15.109, Class B

Approved for and on behalf of Integrated Technologies Ltd by:

Andrew Goodenough



Date 12-JUNE-18

Quality Assurance and Regulatory Affairs Manager