

PRODUCT DATA SHEET

PIG SERUM, STERILE, TRACE-HEMOLYZED

Product Codes:

39526-5

Pack size:

500 mL

Description:

Whole blood is collected from healthy pigs (mixed breed, either sex) at the time of slaughter, refrigerated within 6 hours, processed within 24 hours and frozen at -20 °C. The serum is thawed, pooled, passed through a 0.22 micron filter, bottled and returned to -20 °C.

Physical State:

Frozen liquid

Testing:Hemolysis: $\geq 20 \leq 50$ mg/dL

Total Protein: 7.0-10.0 mg /dL

pH: 7.0-8.5

Sterility: Passed, no growth at 14 days

Packaging, shipping/storage:**Packaging:**

PETG bottle

Storage Temperature:

-20 °C or below

Shipping Conditions:

Dry ice

Expiration:

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

Application Notes:

Reportedly useful in IHC assays and protein purification.

P/N: 74036 Rev 01