

NFIB Polyclonal Antibody

Catalog No.	A15294	Category	Polyclonal Antibodies
Applications	WB	Observed MW	47kDa
Cross-reactivity	Human	Calculated MW	22kDa/33kDa/47kDa/53kDa/55kDa

Immunogen Information

Immunogen	A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human NFIB (NP_005587.2).
Gene ID	4781
Swiss prot	O00712
Synonyms	NFIB; CTF; HMGIC/NFIB; NF-I/B; NF1-B; NFI-B; NFI-RED; NFIB2; NFIB3; nuclear factor I B


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - NFIB Polyclonal Antibody (A15294)