


## Anti-VPS26A Antibody (A82694)

### Specifications:

Name:	Anti-VPS26A Antibody
Description:	Goat polyclonal antibody to VPS26A.
Applications:	ELISA, WB
Reactivity:	Human
Immunogen:	Synthetic peptide corresponding to Human VPS26A (C terminal).
Sequence:	C-SPESQASAEQPEM
Host:	Goat
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Concentration:	100 µg at 0.5 mg/ml.
Product Form:	Liquid
Formulation:	Supplied in Tris Buffered Saline, pH 7.30, with 0.02% Sodium Azide and 0.5% BSA.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.


### Images:



Anti-VPS26A Antibody (A82694) staining (1µg/ml) of Human Kidney lysate (RIPA buffer, 35µg total protein per lane). Primary incubated for 1 hour. Detected by chemiluminescence.

## Anti-VPS26A Antibody (A82694)

Images continued:



HEK293 overexpressing VPS26A (RC205353) and probed with Anti-VPS26A Antibody (A82694) (mock transfection in first lane), tested by Origene.