

Mouse anti RFP-Tag mAb

Catalog No.	AE020	Category	Tag Antibodies
Applications	WB, IHC, IP	Observed MW	27kDa
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	Recombinant protein of human RFP-Tag
Gene ID	
Swiss prot	
Synonyms	

Product information

Source	Mouse
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

Recommended Dilutions

WB 1:2000 -
1:5000

IHC 1:50 -
1:200

IP 1:50 -
1:100



Western blot - Mouse anti RFP-Tag mAb (AE020)