

## Rabbit anti V5-Tag pAb

<b>Catalog No.</b>	AE022	<b>Category</b>	Tag Antibodies
<b>Applications</b>	WB	<b>Observed MW</b>	Refer to Figures
<b>Cross-reactivity</b>		<b>Calculated MW</b>	

### Immunogen Information

<b>Immunogen</b>	A synthetic peptide of V5 epitope
<b>Gene ID</b>	
<b>Swiss prot</b>	
<b>Synonyms</b>	


### Product information

<b>Source</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

### Background

### Recommended Dilutions

WB 1:1000 -  
1:2000



Western blot - Rabbit anti V5-Tag pAb (AE022)