

Eblasakimab Biosimilar - Research Grade

More Information

Concentration: 5mg/ml

isotype: IgG4 - kappa

Host: CHO Cells

Target: IL13RA1

Species Reactivity: Human

Specificity: Detects human IL13RA1

Purification Method: This monoclonal antibody was purified using Protein A

Formulation: Sterile, preservative-free, solution in PBS. BSA and Azide free.

Purity: >95% by SDS-PAGE

Endotoxin: ² 1.0 EU/mg as determined by the LAL method

Storage: This biosimilar is stable when stored at 2-8°C. For longer term storage (> 3 months) it is recommended to store this antibody at -20°C or -70°C.

Applications: Functional Assays

Application Notes: This biosimilar is for research use only (RUO); it is not for diagnostic or therapeutic procedures and cannot be purchased by patients.

Use: Products are for research use only. Not for use in diagnostic or therapeutic procedures.

Delivery Time: 4 weeks