

DCAF11 Polyclonal Antibody

Catalog No.	A15519	Category	Polyclonal Antibodies
Applications	WB	Observed MW	61kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	50kDa/58kDa/61kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human DCAF11 (NP_079506.3).
Gene ID	80344
Swiss prot	Q8TEB1
Synonyms	DCAF11; GL014; PRO2389; WDR23; DDB1 and CUL4 associated factor 11

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

This gene encodes a WD repeat-containing protein that interacts with the COP9 signalosome, a macromolecular complex that interacts with cullin-RING E3 ligases and regulates their activity by hydrolyzing cullin-Nedd8 conjugates. Multiple alternatively spliced transcript variants have been found for this gene.

Recommended Dilutions

WB 1:200 -
1:2000



Western blot - DCAF11 Polyclonal Antibody (A15519)