

HLX Polyclonal Antibody

Catalog No.	A7436	Category	Polyclonal Antibodies
Applications	WB	Observed MW	51kDa
Cross-reactivity	Human	Calculated MW	50kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-280 of human HLX (NP_068777.1).
Gene ID	3142
Swiss prot	Q14774
Synonyms	HLX; HB24; HLX1; H2.0 like homeobox


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - HLX Polyclonal Antibody (A7436)