

Part Number(s): **SF-TS-600**

Description: **Adhesive sealing film for PCR**

Specifications:

Dimensions	<p>TOP VIEW</p>
Tolerances	Overall Length X Width: +/- 0.04"
Corner radii	1/16"

Composition

	Film	Adhesive
Color	Clear	Clear
Material	Polyester	Acrylic
Thickness	51 µm	38 - 46 µm
Operational temperature	-20°C to +120°C	

Quality	
DNase free	Yes
RNase free	Yes
Non-pyrogenic	Yes

Storage
Protect from direct sunlight and UV light
Store dry at 21°C
Designed for single use in PCR thermocyclers

Bovine Spongiform Encephalopathy (BSE) and Transmissible Spongiform Encephalopathy (TSE):

None of these products or raw materials are manufactured from or come into contact with animal material.