

CUL4A Polyclonal Antibody

Catalog No.	A6199	Category	Polyclonal Antibodies
Applications	WB	Observed MW	88kDa
Cross-reactivity	Human, Mouse	Calculated MW	76kDa/87kDa

Immunogen Information

Immunogen	A synthetic peptide of human CUL4A
Gene ID	8451
Swiss prot	Q13619
Synonyms	CUL4A; cullin-4A


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:1000



Western blot - CUL4A Polyclonal Antibody (A6199)