

## CGN Polyclonal Antibody

<b>Catalog No.</b>	A15489	<b>Category</b>	Polyclonal Antibodies
<b>Applications</b>	WB	<b>Observed MW</b>	160kDa
<b>Cross-reactivity</b>	Human	<b>Calculated MW</b>	87kDa/136kDa

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein containing a sequence corresponding to amino acids 904-1203 of human CGN (NP_065821.1).
<b>Gene ID</b>	57530
<b>Swiss prot</b>	Q9P2M7
<b>Synonyms</b>	CGN; cingulin


### Product information

<b>Source</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification method</b>	Affinity purification
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

### Background

### Recommended Dilutions

WB 1:200 -  
1:2000



Western blot - CGN Polyclonal Antibody (A15489)