

HRP-conjugated Mouse anti Myc-Tag mAb

Catalog No.	AE026	Category	Tag Antibodies
Applications	WB	Observed MW	Refer to figures
Cross-reactivity		Calculated MW	

Immunogen Information

Immunogen	A synthetic Peptide of human Myc-Tag
Gene ID	
Swiss prot	
Synonyms	


Product information

Source	Mouse
Isotype	
Purification method	
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:2000 -
1:5000



Western blot - HRP-conjugated Mouse anti Myc-Tag mAb (AE026)