

IGFBP5 Polyclonal Antibody

Catalog No.	A1720	Category	Polyclonal Antibodies
Applications	WB	Observed MW	34kDa
Cross-reactivity	Human, Mouse	Calculated MW	30kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-272 of human IGFBP5 (NP_000590.1).
Gene ID	3488
Swiss prot	P24593
Synonyms	IGFBP5; IBP5; insulin-like growth factor-binding protein 5


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - IGFBP5 Polyclonal Antibody (A1720)