

## Phospho-Tau Thr181 (p-tau181) Polyclonal Antibody

---

### Specification Sheet

Catalog Number	Unit Size
TBS10033-0.1MG	0.1mg
TBS10033-0.5MG	0.5mg

**Description:** Phospho-Tau Thr181 (p-tau181) Polyclonal Antibody. This antibody is validated as capture in ELISA.

**Specificity:** This antibody reacts specifically with phosphor-tau Thr181.

**Host Animal:** Rabbit.

**Isotype:** IgG.

**Immunogen:** Synthesized phosphor-tau thr181 peptide

**Format:** Purified, liquid.

**Purification:** DEAE chromatography.

**Concentration:** 1.0 mg/ml.

**Affinity Constant:** Not determined.

**Buffer:** 10 mM PBS, 0.05% NaN<sub>3</sub>, pH7.2.

**Purity:** >95% by HPLC and SDS-PAGE.

**Application:** For human phosphor-tau Thr181 quantitative assays by EIA or CLIA and Lateral Flow.

**Storage:** Store at -20°C for 2years.

**For Research Use only**