

STARD3 Polyclonal Antibody

Catalog No.	A7776	Category	Polyclonal Antibodies
Applications	WB, IF	Observed MW	45kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	48kDa/49kDa/50kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-445 of human STARD3 (NP_006795.3).
Gene ID	10948
Swiss prot	Q14849
Synonyms	STARD3; CAB1; MLN64; es64; stAR-related lipid transfer protein 3

Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Background

This gene encodes a member of a subfamily of lipid trafficking proteins that are characterized by a C-terminal steroidogenic acute regulatory domain and an N-terminal metastatic lymph node 64 domain. The encoded protein localizes to the membranes of late endosomes and may be involved in exporting cholesterol. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 -
1:2000

IF 1:50 -
1:200



Western blot - STARD3 Polyclonal Antibody (A7776)