

**PRODUCT DATA SHEET****CHICKEN PLASMA, NON-STERILE, SODIUM EDTA****Product Codes**

33141-1

**Pack size:**

500 mL

**Description:**

Healthy animals 6-8 weeks of age, mixed gender and breed; collected on wet ice following USDA guidelines in 2.5% Sodium EDTA solution at 1:10 dilution.

**Physical State:**

Frozen liquid

**Packaging, shipping/storage:****Packaging**

Polypropylene bottles

**Storage Temperature**

-20 °C or below

**Shipping Conditions**

Dry ice

**Expiration**

Product quality is guaranteed to meet Pel-Freez Biologicals' specifications for 1 year from the date of receipt by the customer as long as the product is stored in accordance with the indicated storage conditions.

**Application Notes:**

Reportedly useful for protein purification

P/N: 74042 Rev 01