

PDP2 Polyclonal Antibody

Catalog No.	A14416	Category	Polyclonal Antibodies
Applications	WB, IHC	Observed MW	57kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	59kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 67-250 of human PDP2 (NP_065837.1).
Gene ID	57546
Swiss prot	Q9P2J9
Synonyms	PDP2; PDPC 2; PPM2B; PPM2C2; pyruvate dehydrogenase phosphatase catalytic subunit 2

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background


Recommended Dilutions

WB 1:500 -
1:2000

IHC 1:50 -
1:200



Western blot - PDP2 Polyclonal Antibody (A14416)



Immunohistochemistry - PDP2 Polyclonal Antibody (A14416)