


## Anti-UACA Antibody [UACA/1222] (A249693)

### Specifications:

Name:	Anti-UACA Antibody [UACA/1222]
Description:	Mouse monoclonal [UACA/1222] antibody to UACA.
Specificity:	UACA (Uveal Autoantigen with Coiled-coil domains and Ankyrin repeats) is a 1,416 amino acid nuclear membrane protein. It was originally identified as an autoantigen in patients with panuveitis, a characteristic of Vogt-Koyanagi-Harada disease, and in patients with Graves' disease. UACA was also later identified as Nucling, an mRNA differentially expressed in F9 embryonal carcinoma cells during cardiac muscle differentiation. UACA appears to function as a pro-apoptotic protein that recruits the apaf-1-pro-caspase-9 complex for the induction of apoptosis to mediate the cell-death pathway.
Applications:	WB
Recommended Dilutions:	WB: 1-2 µg/ml
Reactivity:	Human, Mouse
Immunogen:	Recombinant human UACA protein.
Host:	Mouse
Clonality:	Monoclonal
Clone ID:	UACA/1222
Isotype:	IgG1
Light Chains:	kappa
Conjugate:	Unconjugated
Purification:	Protein A/G chromatography.
Concentration:	200 µg/ml
Product Form:	Liquid
Formulation:	Supplied in 10mM Phosphate Buffered Saline with 0.05% BSA and 0.05% Sodium Azide.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
General Notes:	This monoclonal antibody is also available in a different formulation without BSA and Sodium Azide - Anti-UACA Antibody [UACA/1222] - BSA and Azide free (A252873).
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

## Anti-UACA Antibody [UACA/1222] (A249693)

Images:



Western blot analysis of A549 cell lysate using Anti-UACA Antibody [UACA/1222].