

NAA20 Polyclonal Antibody

Catalog No.	A7105	Category	Polyclonal Antibodies
Applications	WB, IHC, IF, IP	Observed MW	20kDa
Cross-reactivity	Human, Mouse, Rat	Calculated MW	12kDa/20kDa


Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-178 of human NAA20 (NP_057184.1).
Gene ID	51126
Swiss prot	P61599
Synonyms	NAA20; NAT3; NAT3P; NAT5; NAT5P; dj1002M8.1; N-alpha-acetyltransferase 20


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background



Immunofluorescence - NAA20 Polyclonal Antibody (A7105)



Immunoprecipitation - NAA20 Polyclonal Antibody (A7105)


Recommended Dilutions

WB 1:500 -
1:2000


IHC 1:50 -
1:200

IF 1:50 -
1:100

IP 1:50 -
1:200



Western blot - NAA20 Polyclonal Antibody (A7105)



Immunohistochemistry - NAA20 Polyclonal Antibody (A7105)