

AMBRA1 Polyclonal Antibody

Catalog No.	A12578	Category	Polyclonal Antibodies
Applications	WB	Observed MW	175kDa
Cross-reactivity	Mouse	Calculated MW	84kDa/132kDa/135kDa/139kDa/142kDa

Immunogen Information

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 961-1208 of human AMBRA1 (NP_060219.2).
Gene ID	55626
Swiss prot	Q9C0C7
Synonyms	AMBRA1; DCAF3; WDR94; autophagy and beclin 1 regulator 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - AMBRA1 Polyclonal Antibody (A12578)