



DATA SHEET

HUMAN PNEUMADIN

Catalog no.:	AE1134AG.1 / AE1134AG.2
Immunogen:	Human Pneumadin (AGEPKLDAGV)
Reference no.:	P22103
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0,1% acetic acid for required concentration
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Only for research use

N:\WINWORD\DAT\DATEN\LAETTER\ANTIKOERPER\MISCELLANEOUS\AE 1134.DOC