

IGFBPL1 Polyclonal Antibody

Catalog No.	A16635	Category	Polyclonal Antibodies
Applications	WB	Observed MW	35kDa
Cross-reactivity	Mouse, Rat	Calculated MW	29kDa

Immunogen Information

Immunogen	A synthetic peptide of human IGFBPL1
Gene ID	347252
Swiss prot	Q8WX77
Synonyms	IGFBPL1; IGFBP-RP4; IGFBPRP4; bA113O24.1; insulin-like growth factor-binding protein-like 1


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Background

Recommended Dilutions

WB 1:500 -
1:2000



Western blot - IGFBPL1 Polyclonal Antibody (A16635)