

Neurofilament L (NFL) Mouse Monoclonal Antibody

Specification Sheet

Catalog Number	Unit Size
TBS10027H-0.1MG	0.1mg
TBS10027H-0.5MG	0.5mg

Description: Neurofilament Light (NFL) Mouse Monoclonal Antibody conjugated with HRP. This antibody as a detection antibody can be paired with TBS10025 or TBS10028 as capture.

Specificity: This antibody reacts with NFL protein.

Host Animal: Balb/c mouse.

Isotype: IgG1.

Immunogen: Recombinant NFL protein.

Format: Purified, liquid.

Conjugate: HRP

Concentration: 1.0 mg/ml.

Affinity Constant: Not determined.

Buffer: 10 mM PBS, 1%BSA, 50% glycerol, pH7.2.

Pairing Assays: TBS10025 and TBS10028.

Application: For human NFL quantitative assays by EIA or CLIA and Lateral Flow.

Storage: 4°C for 1 year.

For Research Use only