

## Anti-SYK Antibody (A96427)

### Specifications:

Name:	Anti-SYK Antibody
Description:	Rabbit polyclonal antibody to SYK.
Applications:	WB, IHC, ELISA
Recommended Dilutions:	WB: 1:500-1:1000, IHC: 1:50-1:100, ELISA: 1:5000
Reactivity:	Human, Mouse, Rat
Immunogen:	Synthetic peptide derived from human SYK (amino acids 491-540).
Host:	Rabbit
Clonality:	Polyclonal
Isotype:	IgG
Conjugate:	Unconjugated
Purification:	Purified from rabbit serum by antigen affinity chromatography using the immunizing peptide.
Molecular Weight:	72kDa
Product Form:	Liquid
Formulation:	Supplied in Phosphate Buffered Saline (without Mg2+ and Ca2+), pH 7.4, with 150mM NaCl, 0.02% Sodium Azide, and 50% Glycerol.
Storage:	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
Disclaimer:	This product is for research use only. It is not intended for diagnostic or therapeutic use.

### Images:



Western blot analysis of lysates from COS using Anti-SYK Antibody. The right hand lane represents a negative control, where the antibody is blocked by the immunising peptide.

# antibodies

## Anti-SYK Antibody (A96427)

#### Images continued:



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using Anti-SYK Antibody. The right hand panel represents a negative control, where the antibody was pre-incubated with the immunising peptide.



Western blot analysis of various cells using Anti-SYK Antibody.



Immunohistochemical analysis of paraffin-embedded human breast cancer using Anti-SYK Antibody 1:100 (4°C overnight). The right hand panel represents a negative control, where the antibody was pre-incubated with the immunising peptide.