

ZNF217 Polyclonal Antibody

Catalog No.	A11812	Category	Polyclonal Antibodies
Applications	WB	Observed MW	120kDa
Cross-reactivity	Human	Calculated MW	115kDa

Immunogen Information

Immunogen	A synthetic peptide of human ZNF217
Gene ID	7764
Swiss prot	O75362
Synonyms	ZNF217; ZABC1; zinc finger protein 217


Product information

Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

Recommended Dilutions

WB 1:200 -
1:500



Western blot - ZNF217 Polyclonal Antibody (A11812)