antibodies

Anti-PAG1 Antibody [PAG-C1] (A85668)

Specifications:

Anti-PAG1 Antibody [PAG-C1]
Mouse monoclonal [PAG-C1] antibody to PAG1.
The antibody PAG-C1 recognizes an epitope located in the C-terminal domain of Csk-binding protein (Cbp / PAG), a 46 kDa ubiquitously expressed transmembrane adaptor protein present in membrane rafts (glycosphingolipid-enriched microdomains), which however migrates on SDS PAGE gels anomalously as an 80 kDa molecule.
IP, WB, IHC-P
Mouse, Human, Rat, Bovine
C-terminal peptide (last 15 amino acids) of human Csk binding protein coupled to KLH.
Mouse
Monoclonal
PAG-C1
lgG2b
Unconjugated
Protein A chromatography.
1 mg/ml
> 95% (by SDS-PAGE).
Liquid
Supplied in Phosphate Buffered Saline, pH 7.4, with 15 mM Sodium Azide.
Shipped at 4° C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
This product is for research use only. It is not intended for diagnostic or therapeutic use.

antibodies

Anti-PAG1 Antibody [PAG-C1] (A85668)

Images:



Immunohistochemistry staining (paraffin sections) of PAG/Cbp in germinal center of lymphatic follicle and in dispersed T cells in appendix tissue by Anti-PAG Antibody (A85668). Positive signal in T cells.



Immunoprecipitation of human PAG/Cbp from the lysate of RAJI human Burkitt lymphoma cell line. Western blot was immunostained with Anti-PAG Antibody (A86214). Note: PAG/Cbp is a 46 kDa adaptor protein, which however migrates on SDS PAGE gels anomalously as an 80 kDa molecule. Lane 1 & 2: Immunoprecipitation with Anti-PAG Antibody (A85668). Lane 3: Immunoprecipitation with Anti-PAG Antibody (polyclonal antibody).