

CCIN Polyclonal Antibody

Catalog No.	A3777	Category	Polyclonal Antibodies
Applications	WB	Observed MW	72kDa
Cross-reactivity	Human	Calculated MW	66kDa

Immunogen Information

Immunogen	A synthetic peptide of human CCIN
Gene ID	881
Swiss prot	Q13939
Synonyms	CCIN; BTBD20; KBTBD14; calicin

Product information


Source	Rabbit
Isotype	IgG
Purification method	Affinity purification
Storage	Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, pH7.3.

Background

The protein encoded by this gene is a basic protein of the sperm head cytoskeleton. This protein contains kelch repeats and a BTB/POZ domain and is necessary for normal morphology during sperm differentiation. This gene is intronless.

Recommended Dilutions

WB 1:500 -
1:1000



Western blot - CCIN Polyclonal Antibody (A3777)